(M)ELISSA

16780011

DK	Lysterapilampe / LYSTERAPIENHED	2
SE	Dagsljuslampa/ljusterapienhet	11
NO	Dagslys Lampe / LYSTERAPIENHET	20
FI	Päivänvalolamppu / VALOHOITOYKSIKKÖ	29
UK	Daylight lamp / PHOTOTHERAPY UNIT	38
DE	Tageslichtlampe FOTOTHERAPIEGERÄT	47

INDLEDNING

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INTRODUCTION

To get the best out of your new appliance, please read this user guide carefully before using it for the first time. Take particular note of the safety precautions. We also recommend that you keep the instructions for future reference, so that you can remind yourself of the functions of your appliance.

Expectation

This phototherapy unit is intended for light therapy to treat the seasonal affective disorder (S.A.D.) It is a convenient and effective way of compensating for the lack of light without resorting to medication in autumn and winter.

SAFETY MEASURES

General

- Incorrect use of this appliance may cause personal injury and damage.
- Use for its intended purpose only. The manufacturer is not responsible for any injury or damage resulting from incorrect use or handling (see also Guarantee Terms).
- For domestic use only. Not for outdoor or commercial use.
- For connection to 230 V, 50 Hz only.
- Do not immerse the appliance, cord or plug in water and make sure no water enters the appliance.
- Do not use in humid surroundings.
- Never touch the appliance, cord or plug with wet or damp hands.
- Do not expose to direct sunlight, high temperatures, humidity, noise or corrosive substances.

- Keep the appliance under constant supervision while in use. Keep an eye on children.
- Do not cover.
- Use under temperature 10-25 °C.
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Children should be supervised to ensure that they do not play with the appliance.

Cord and plug

- Keep the cord and appliance away from heat sources, hot objects and naked flames.
- Check that it is not possible to pull or trip over the cord or any extension cord.
- Switch off the appliance; remove the plug from the socket when the appliance is not in use and before cleaning or replacement of the fluorescent tubes. Avoid pulling the cord when removing the plug from the socket. Instead, hold the plug.
- Handle fluorescent tubes carefully as they may splinter and cause injuries to your hands.
- Ensure that the cord is fully extended.
- The cord must not be twisted or wound around the appliance.


- Check regularly for damage to appliance, cord and plug. Do not use if any parts are damaged or if the appliance has been dropped or damaged in any way whatsoever.
- If the appliance, cord or plug has been damaged, have them inspected and if necessary repaired by an authorized repair engineer. Never try to repair the appliance yourself. Please contact the store where you bought the appliance for repairs under guarantee. Unauthorized repairs or modifications to the appliance will invalidate the guarantee.

Using the phototherapy unit

- Do not use in rooms without some other light source. For use as a light supplement only, not illumination.
- Do not look directly at the appliance for long periods. Follow the user guide below.
- Initial use may cause tired eyes and/or headache. This is normal, and such effects will soon disappear. If symptoms persist, consult your doctor.
- Consult your doctor before using this appliance if you
 - o suffer or have suffered from severe depression
 - o are hypersensitive to light
 - o suffer from any type of eye disease
 - o use particular forms of medicine such as psychopharmacia or malaria pills
 - o suffer from high blood pressure.
- Consult your doctor, if the use of the appliance results in hyperactivity, exhaustion, depression, restlessness, difficulties in sleeping and/or other adverse effects.

KEY TO MAIN COMPONENTS OF THE APPLIANCE

1. Screen
2. Brackets
3. Mounting hole
4. Carrying handle
5. Fold-out stand
6. On/off switch
7. Retaining clips
8. Reflector
9. Tubes, 2 x 36W

USING THE APPLIANCE

The appliance is used for light therapies that counteract winter fatigue, which can cause the following symptoms:

- fatigue
- bad mood or moodiness
- lack of energy
- increased need for sleep
- lack of concentration
- lack of appetite

Over 60% of the population is believed to experience one or more symptoms of winter blues due to the low intensity of the winter sun.

The most effective time of day to take a light therapy is between 6 am and 8 pm for a recommended period of 2 hours per day.

You can take light therapies as often as you like but the treatment is most effective if you take light therapies of the below-mentioned duration at least five days in a row.

Positioning the appliance

- Unfold the stand at the back of the appliance and place the appliance on a table.

or

- Mount the appliance on the wall using the enclosed wall brackets. Check before mounting that the wall and the screws can carry the weight of the appliance.

Use

- Insert the plug in the mains socket and switch on the power.
- Switch on the on/off button.
- Sit down in front of the appliance in a position where you can comfortably see the light and with your face 20-60 cm from the screen. You do not need to sit directly in front of the screen.
- The duration of the light therapy depends on your distance from the screen:

Distance to screen	Max duration of light therapy
60 cm	Approx. 2 hours
40 cm	Approx. 1 hour
20 cm	Approx. ½ hour

- The light therapy can be interrupted and you can do other things at the same time (e.g. read or watch television) as long as you sit within 60 cm from the appliance.
- Do not look directly into the light constantly, but look at it regularly.


- Always turn off the appliance using the on/off button, and unplug it, when it is not in use.
- Allow the appliance to cool for at least 15 minutes before putting it away or moving it.

TIP: Carry the appliance by the carrying handle on the top of the rear part.

REPLACEMENT OF FLUORESCENT TUBES AND STARTING SWITCH

The lifetime of the fluorescent tubes is approx. 8,000 hours of normal use. If a tube fails, it must be replaced.

- **IMPORTANT!** Discard used fluorescent tubes in accordance with the rules for environmentally harmful waste, e.g. by taking them to a waste recovery site.
- **IMPORTANT!** Always make sure to switch off the appliance and remove the plug from the mains socket before taking it apart to avoid the risk of electric shocks.
- **IMPORTANT!** Always let the appliance cool down before changing the tubes.
- Place the appliance on a non-scratching surface, screen down.

<p>1. Slacken the retaining screws on the back using a screwdriver, lay the apparatus on its back and remove the screen.</p> <p> Note</p> <p>The two brackets on the front of the screen will become loose when the screws are slackened. Take care not to lose them, or the screen cannot be refitted correctly.</p>	
<p>2. Bend the tube's four transparent retaining clips by pulling them gently upwards.</p> <p> Note</p> <p>Cable ties are used for transport protection only.</p>	
<p>3. Carefully pull the tube up over the housing rim (1) and out of the housing, then sideways out of the frame.</p>	
<p>4. When inserting the tube, following these instructions in reverse order.</p>	

CLEANING

When cleaning the appliance, you should pay attention to the following points:

- Remove plug from the mains socket, and allow the appliance to cool down before cleaning it.
- Never submerge the apparatus in water, and no water must be allowed to get into it.

External cleaning

- Clean the appliance by wiping it with a damp cloth. A little detergent can be added if the appliance is heavily soiled.
- Never use a scouring pad, steel wool or any form of strong solvents or abrasive cleaning agents to clean the appliance, as those may damage the outside surfaces of the appliance.
- **IMPORTANT!** Allow to dry fully after cleaning before using again. Do not use appliance if damp.

Internal cleaning

- Remove the screen as described above, wipe off reflector, inside of screen and tubes with a dry, soft cloth. Do not use any kind of cleaning agents.

SPECIFICATIONS

Voltage	220–240V
Wattage	72W
Hz	50
LUX	10,000 at 15cm
Lamp type	PLL 36W 2G11 4pin H
Fuse	T3.15AH250V


Operating conditions:

Temperature: +10 °C to 25 °C
Relative Humidity: 30 % to 90 %

Storage conditions:

Temperature: +5 °C to 50 °C
Relative Humidity: 30 % to 90 % (no condensation)
Class I appliance; IPX0

EXPLANATION OF SYMBOLS ON THE EQUIPMENT

	Consult instructions for use
	Caution
	If you want to dispose this product, do not mix with general household waster. There is a separate collection system for used electronics products in accordance with legislation under WEEE Directive (Directive 2002/96/EC) and is effective only with in European Union.
	Manufacturer
	Date of manufacture
	Serial number
	This appliance meets the requirement of Medical Device Directive 93/42/EEC for medical appliances.
	“ON” / “OFF”
	Authorized Representative in the European Community

EMC Information-Guidance and Manufacture's Declaration
Technical Description Concerning Electromagnetic Emission

Table 1: Declaration - electromagnetic emissions

Guidance and manufacturer's declaration – electromagnetic emissions

The “Daylight” is intended for use in the electromagnetic environment specified below. The customer or the user of the “Daylight” should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The “Daylight” uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The “Daylight” is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000–3–2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000–3–3	Complies	

Table 2: Guidance & Declaration – electromagnetic immunity

Guidance and manufacturer's declaration – electromagnetic immunity

The “Daylight” is intended for use in the electromagnetic environment specified below. The customer or the user of the “Daylight” should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000–4–2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000–4–4	±2 kV for power supply lines ±1 kV for input/ output lines	±2 kV for power supply lines ±1 kV for input/ output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000–4–5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000–4–11	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the “Daylight” requires continued operation during power mains interruptions, it is recommended that the “Daylight” be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000–4–8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE UT is the a.c. mains voltage prior to application of the test level.

Table 3: Guidance & Declaration – electromagnetic immunity concerning Conducted RF & Radiated RF

Guidance and manufacturer's declaration – electromagnetic immunity			
The “Daylight” is intended for use in the electromagnetic environment specified below. The customer or the user of the “Daylight” should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000–4–6	3 Vrms 150 kHz to 80 MHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of the “Daylight”, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=1.2\sqrt{P}$ $d=1.2\sqrt{P}$ 80MHz to 800MHz $d=2.3\sqrt{P}$ 800MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000–4–3	3 V/m 80 MHz to 2,5 GHz	3 V/m	
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the “Daylight” is used exceeds the applicable RF compliance level above, the Medical Daylight should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the “Daylight”.			
b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.			

Table 4: Recommended separation distances between portable and mobile RF communications equipment and the model

Recommended separation distances between portable and mobile RF communications equipment and the Medical Daylight			
The “Daylight” is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Medical Daylight can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the “Daylight” as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = [\frac{3,5}{V_1}] \sqrt{P}$	80 MHz to 800 MHz $d = [\frac{3,5}{E_1}] \sqrt{P}$	800 MHz to 2,5 GHz $d = [\frac{7}{E_1}] \sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Electromagnetic fields (EMF)

This appliance complies to all standards regarding electromagnetic fields (EMF). If handled properly and in accordance with the instructions given in this manual, the appliance is safe to use, based on today’s scientific knowledge.

FOR UNITED KINGDOM ONLY

Plug wiring:

This product is fitted with a BS 1363 13-amp plug. If you have to replace the fuse, only those that are ASTA or BSI approved to BS1362 and with a rated current of 13 amps should be used. If there is a fuse cover fitted, this cover must be refitted after changing the fuse. If the fuse cover is lost or damaged, the plug must not be used. You must also check if the socket outlets in your home fit with the plug of the appliance. If the socket outlet in your home does not fit with the plug, the plug must be removed and disposed of safely as insertion of the plug into the socket is likely to cause electric hazard.

As this product is a class II appliance, the pin of BS plug for earth connection is false.

INFORMATION ON DISPOSAL AND RECYCLING OF THIS PRODUCT

Please note that this Adexi product is marked with this symbol: 

This means that this product must not be disposed of along with ordinary household waste, as electrical and electronic waste must be disposed of separately.

In accordance with the WEEE directive, every member state must ensure correct collection, recovery, handling and recycling of electrical and electronic waste. Private households in the EU can take used equipment to special recycling stations free of charge. In certain member states, used apparatus can be returned to the dealer where they were bought on the condition you buy new products. Contact your retailer, distributor or the municipal authorities for further information on what to do with electrical and electronic waste.

GUARANTEE TERMS

The guarantee does not apply:

- if the above instructions are not followed
- if the appliance has been interfered with
- if the appliance has been mishandled, subjected to rough treatment, or has suffered any other form of damage
- if faults have arisen because of faults in your electricity supply.

Due to the constant development of our products in terms of function and design, we reserve the right to make changes to the product without prior warning.

FREQUENTLY ASKED QUESTIONS

If you have any questions regarding use of the appliance and cannot find the answer in this user guide, please try our website at www.adexi.eu.

Go to the "Consumer Services" menu, click on "Frequently Asked Questions" to view the most frequently asked questions.

You can also see contact details if you need to contact us for technical issues, repairs, accessories and spare parts.

IMPORTER

Adexi Group
www.adexi.eu

Manufacturer

Zero-Plus International Ltd.
Room 1004, 10/F., Join-In Hang Sing Centre,
71-75 Container Port Road, Kwai Chung,
New Territories, Hong Kong

European (EU) Authorized Representative

Advena Ltd.
Thome Widgery House,
33 Bridge Street,
Hereford HR4 9DQ
UK

We cannot be held responsible for any printing errors.

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